



Your Guide to

NEOCAREPAL

An all-in-one calculator app for NICU care



NeoCarePal is your trusted calculator companion, providing a single destination for:







- Helpful Tip -


With a fully functional "offline" mode, the app's calculators do not require an internet connection for use.


NICU=neonatal intensive care unit.

How NeoCarePal Can Help

Do you find yourself juggling multiple websites, references, and clinical tools as you care for your patients? NeoCarePal can help make this simpler by providing multiple calculators in one place.







Drawing from established guidelines and clearly referenced scientific evidence, NeoCarePal allows access to clinical calculators with the convenience to:



Input patient metrics and parameters



Tab between different calculators applying the same patient metrics



Save anonymized patient data for up to 24 hours



Use built-in calculator tools in the offline mode



Streamline calculator inputs with up to 10 customizable preset profiles



Access calculators and reference materials at third-party websites through the search tool

NeoCarePal FAQs

How do I download the app?

For iPhone users, the NeoCarePal app is available for download from the App Store. For Android users, NeoCarePal is available for download on Google Play.

Do I need to create an account?

Yes, to register as a user of NeoCarePal, you will have to input your first and last name, your healthcare role, and the hospital name. You will be prompted to enter your email and create a password, and then use your email and password each time you log in.

How does the NeoCarePal app manage identifiable patient information?

All identifiable patient information is stored locally on your device. No identifiable patient information is collected, transmitted, used, or accessed in any way by Chiesi USA. As a result, while the NeoCarePal app's security conforms with industry standards, you are responsible for any identifiable patient information entered into the app. The information you enter into the NeoCarePal app is not integrated into the patient's electronic health record. When you create a customized preset profile, you can input patient metrics that can be locally stored on the app. NeoCarePal generates a unique NCP code for that profile. You can use the NCP code to connect the data to the patient.

Chiesi USA does not have a HIPAA "business associate" arrangement with your hospital regarding the NeoCarePal app, nor is Chiesi USA able to review your hospital's HIPAA policies on use of apps of this nature. Rather, we are making the NeoCarePal app available directly to you.

How long can I store data in the NeoCarePal app?

Information can be stored in the app for up to 24 hours, allowing you to track inputs over that 24-hour period. Twenty-four hours after each calculation is performed, the data will be erased from the app.

Download NeoCarePal Today



Practical Integration

NeoCarePal seamlessly integrates various neonatal care calculators and clinical practice tools into a single, user-friendly app.



Streamlined Workflow

By reducing the need to juggle multiple references and resources, NeoCarePal may simplify the day-to-day work of NICU staff, saving time and reducing cognitive load.



Patient Care Support

With NeoCarePal, healthcare providers can instantly access and rely on calculations and relevant information, providing dependable support that is essential to quality neonatal care.

Download the NeoCarePal App

iPhone Users

Download at bit.ly/NeoCarePalApp

Android Users

Download at bit.ly/NeoCarePalGoogle

